






MAXX-1100

CT ROOM LASER

By Cemar Electro Inc.
www.cemarelectro.com



— CEMARLINE —



CEMARLINE



MAX-1100

FEATURING:

- WINDOWS BASED
- DICOM 3.0 RT
- ETHERNET AND ZIGBEE WIRELESS COMMUNICATION
- SYMBOL BASED REMOTE CONTROL


- ACCURACY +/- 0.2mm
- RANGE OF TRAVEL +/- 300mm

WAVELENGTH:

RED 635nm

GREEN 520nm

BLUE 450nm



Canada: 

Cemar Electro Inc.
1370 55th Ave.
Lachine, QC
Canada H8T 3J8

Tel: 514-631-5807
Fax: 514-631-7505
Toll Free: 1-800-298-5273
cemar@cemarelectro.com

United States: 

Cemar Electro Inc.
100 Walnut St.
Champlain, NY
USA 12919

Tel: 518-298-3065
Fax: 518-298-2927
Toll Free: 1-800-298-5273
cemarusa@cemarelectro.com

