

LIBERTY DERMOGRAPH

Precise layout for targeted marking of irradiation fields

**CORD
LESS**

- ✓ Simple, intuitive handling
- ✓ A lightweight, ergonomic pen
- ✓ Cordless pen for complete freedom
- ✓ Pain-free implantation for the patient
- ✓ Sterile, single-use, retractable needles for patient and operator safety

- ✓ **COMFORTABLE**
- ✓ **FINE**
- ✓ **LIGHTWEIGHT**



Liberty's plus points

The Liberty dermograph is a lightweight, ergonomic dermograph, ideal for marking radiotherapy points.

Thanks to its needles, RadSafe® pigment implantation is painless for the patient and extremely precise.

Its simple, intuitive operation enables operators to quickly get to grips with the Liberty.

Use

Marking of radiotherapy points.

Marking speed: **5 - 5.5 V**

Needle exit for marking: **min. 1 mm**

EASY ADJUSTMENT OF NEEDLE DEPTH

Rotating system



Liberty box components



Technical specifications

Materials	Matt aluminum
Pen size	18.1 cm (with needle)
Pen weight	94 g
Battery weight	42 g
Pen diameter	1.8 cm
Power	5 volt
Battery capacity	900 MA (x2)
Battery life	4 - 5 hr (x2)
Charging time	2 - 3 hr
Other characteristics	<ul style="list-style-type: none"> ● Powerful ● Silent ● Vibration -free

