

# easyBreast

## Breast Board


eaSyBreast is an innovative breast board made of carbon fiber, extremely lightweight, and robust. This device allows to manage the angulation of the main board and the height of the arms separately to prevent collisions and minimize shoulder extension.

The position of the armrests and wrist supports can be managed independently, and the Velcro system contributes to patient comfort. Additionally, they include cushion markers to enhance immobilization cushion fixation and repositioning.

The easyBreast II version also incorporates an anchoring system for indexing breast multipoint masks.



**Tilting System:**  
independent  
adjustment of arm  
and board height.



Product	Reference
<b>Immobilisation System</b>	
eaSyBreast	AG1001048
eaSyBreast II	AG1001072
<b>Accessories</b>	
Slide stopper	AG1001062