

eXaHip

Immobilisation System for the Pelvis

Combines the advantages and functionality of the footrest (eXaFeet) and the knee support (eXaKnee). The combination of eXaFeet and eXaKnee along with the eXaHip base enhances and facilitates patient repositioning.

eXaFeet can be indexed along the eXaHip base in up to 10 different positions (A-J) and can be adjusted at 4 different angles (0°, 5°, 10°, 15°).

eXaKnee can be indexed in 3 different positions on the eXaHip base, which can be indexed to the treatment or CT couch by two indexing bars.

It includes shoulder retractors as an optional accessory, which can be fixed in up to 5 positions (A-E) depending on patient height.

eXaHip is compatible with HexaPOD and it is also available in an MR version.

“eXaHip is designed to improve the precision of positioning and repositioning of the lower abdominal region and the extremities”

Product	Reference
Immobilisation System	
eXaHip	AG1001020-AHEXA
eXaHip MR	AG1001020-AHEXA-MR
Accessories	
eXaHip shoulder retractors	1020.B13

