

# Suremark

## Precision Skin Marking and Labeling

### Suremark labels

Suremark are one of the most popular and versatile adhesive markers. A small tab ensures easy removal. Available in different sizes according to size of the ball.



Reference	Description	Quantity
SL-10	Ball 1 mm of adhesive label 15 mm	110 per box
SL-15	Ball 1.5 mm of adhesive label 15 mm	110 per box
SL-20	Ball 2 mm of adhesive label 15 mm	110 per box

### Suremark Tabs

Suremark tabs are designed to maintain all the characteristics of Suremark Labels with the added benefit of a strip for prominent injuries which adapts to the skin. Available in different sizes, according to the size of the ball.



Reference	Description	Quantity
SL-15T	Ball 1.5 mm of adhesive label 35 mm	110 per box
SL-20T	Ball 2 mm of adhesive label 35 mm	110 per box
SL-25T	Ball 2.5 mm of adhesive label 35 mm	110 per box

### Suremark Wire

The most popular adhesive wire-shaped marker from Suremark. Ideal for marking long injuries or scars. Available in different sizes, according to the wire diameter.



Refer.	Description	Quantity
SL-W3	0.3 mm $\varnothing$ Diamond-shaped adhesive covered wire	336cm per box
SL-W8	0.8 mm $\varnothing$ Diamond-shaped adhesive covered wire	336cm per box
SL-W15	1.5 mm $\varnothing$ Diamond-shaped adhesive covered wire	336cm per box

### Suremark Visionmark


Specially formulated as a non-metallic marker, it allows for precision in various applications. Available in different shapes, according to the ball.



Reference	Description	Quantity
V-20	Ball 2.0 mm of adhesive label	110 per box
V-25	Ball 2.5 mm of adhesive label	110 per box
V-30	Ball 3 mm of adhesive label	50 per box

### Suremark CT Mark

This adhesive marker is manufactured from nonmetallic materials, the result of years of research. The mark generates few artefacts and is perfect for nearly all CT marking applications. Available in different shapes, according to the size of the ball or wire diameter.



Refer.	Description	Quantity
CT-23	Ball 2.3 mm of adhesive label	110 per box
CT-W10	1 mm $\varnothing$ Diamond-shaped adhesive covered wire	330cm per box
CT-W20	2 mm $\varnothing$ Diamond-shaped adhesive covered wire	330cm per box

### Suremark PortalMark Labels

Instead of permanent or temporal tattoos, Suremark have developed PortalMark™ adhesive markers made from breathable, water-resistant materials that stay on the skin up to two weeks without peeling off.



Reference	Description	Quantity
PM-150	Isocentre Markers and field lines	50 sheets per box
PM-325	Adhesive Isocentre Markers	110 per box