

RADSAFE® Ink

RADSAFE® Sterile Tattoo Ink

This product contains permanent ink, sterile, hypoallergenic and non-toxic. The ink cartridge is sealed and radiated with gamma radiation of 25 Gy to guarantee sterility and offer improved quality and patient care. It is class II b medical device.

The cartridge is designed for easy handling and standardisation of tattoo methods. The ink flow is easy to adjust to avoid spilling or dispersing.

The ink does not turn blue with the passage of time, on the contrary, the colour fades. It can be applied directly to the zone or through the cup and lancet provided in the kit.

Main Characteristics:

- ▶ Sterile, non-toxic and hypoallergenic.
- ▶ Ready and easy to use.
- ▶ Available colours: white, black, red, green.
- ▶ Composition available for patients with allergies.
- ▶ Allows for standardisation of the tattoo process.



RadSafe Black
Kit



RadSafe Black
Monodose



RadSafe White
Kit



Available Formats

Producto	Reference	Characteristics	Presentation
Radsafe Black Kit	AG10040121.N	Kit (monodose dispenser, lancet and cup)	50 kits/bag
Radsafe Black Ink	AG10040122.N	Monodose dispenser	100 pcs/bag
Radsafe White Kit	AG10040121.B	Kit (monodose dispenser, lancet and cup)	50 kits/bag
Radsafe White Ink	AG10040122.B	Monodose dispenser	100 pcs/bag