

EXTREMITIES


eaSyFix

Limbs Immobilization Device


eaSyFix

Limbs Immobilization Device

eaSyFix is an immobilization system designed to provide precise immobilization of limbs and other challenging anatomical areas. Equipped with cushion adapters to secure immobilization cushions and multiple indexing areas for thermoplastic masks, eaSyFix ensures a tailored and accurate immobilization for each patient's unique anatomy. This level of customization guarantees unparalleled precision when repositioning during limb treatments.



Versatile, Customizable & Accurate:



Legs



Arms



Back



Chest

Features:

- Compatible with anchor push-pins masks.
- Manufactured of non-artefactant and low attenuation materials.
- Includes 5 indexing cushions adapters.
- Its extensive configuration options make it a high versatile device that can be used not only for immobilizing limbs but also in challenging areas with great precision, such as the chest, breast, or pelvis.
- Fully indexable to any simulation or treatment couch with standard indexing bars.

eXaCast®, High Accuracy Thermoplastic Masks.

Reference	Thermoplastic	Dimensions	Immobilization area
AGEX21123-MC	Microperfo 2,3 mm thickness	150x130 mm	Arm
AGEX21123L-MC	Microperfo 2,3 mm thickness	165x195 mm	Knee
AGEX41123-MC	Microperfo 2,3 mm thickness	355x172,5 mm	Leg
AGEX31123-MC	Microperfo 2,3 mm thickness	262,5x220 mm	Feet
AGP41123K-MC	Microperfo 2,3 mm thickness	300x344 mm	Chest Pelvis
AGB21123-MC	Microperfo 2,3 mm thickness	250x430 mm	Chest Breast
AGBCP1123	Microperfo 2,3 mm thickness	520x450 mm	Chest Abdomen



Moldcare® cushions

Reference	Dimensions
Moldcare R1 II BR-1	40 x 60 cm
Moldcare R1 II BR-3	60 x 60 cm



eXaSkin Plus®

Reference	Products
AG1001003-150	eXaSkin Plus 150
AG1001003-225	eXaSkin Plus 225
AG1001003-650	eXaSkin Plus 650 Breast





Anatomical Geometry, S.L.

Parque Industrial y de Servicios Aljarafe.
C/ Formación N° 6, CP 41927,
Mairena del Aljarafe, Sevilla.

Tel. switchboard: +34 911 436 793.
Fax.: +34 911 238 549.

info@anatge.com

www.anatge.com



FONDO EUROPEO DE
DESARROLLO REGIONAL



Unión Europea

UNA MANERA
DE HACER EUROPA