


eaSyFrame SFV

Immobilisation system for head & shoulders compatible with Multi-point mask with Push-pins



Main Features:

- Width (53 cm) to prevent collisions and compatible with HeXaPOD.
- Ultra-light equipment, made of composite material for minimal attenuation.
- Compatible with standard TIMO (Foam) and SILVERMAN (PVC) headrests.
- Two different mask fixation positions.
- MR Conditional up to 3 Tesla.
- Low-density version available (eaSyFrame S).

Product	Reference
Immobilisation System	
eaSyFrame SFV-MR ultra-light	AG1001052.1
eaSyFrame S	AG1001027

eaSyFrame SLFV

Immobilisation system for head & shoulders compatible with multi-point mask with mush-pins and shoulders retractor optionally



Main Features:

- Width (53 cm) to prevent collisions and compatible with HeXaPOD.
- Lightweight, made of low-density materials for minimal attenuation.
- Compatible with standard TIMO (Foam) and SILVERMAN (PVC) headrests.
- Two different mask fixation positions.
- MR Conditional up to 3 Tesla.
- Shoulder retractor with 10 positions (optional).

Product	Reference
Immobilisation System	
eaSyFrame SLFV-MR	AG1001052
Shoulders retractor	AG1001056