

iCe Arm Support

Thorax immobilisation system MR compatible

iCe Arm Support is a thorax and arms positioning system.

Main Features:

- Handgrips with 10 different positions (1-10), enabling adjustment to the patient's needs.
- Compatible with standard TIMO (Foam) and SILVERMAN (PVC) headrests or Donut.
- Includes two cushion markers for MOLDCARE type cushions.
- indexable to the treatment or simulation couch.
- MR CONDITIONAL up to 3 Tesla



Product	Reference
Immobilisation System	
iCe Arm Support	AG1001074