

WingFrame L

Immobilisation system for head & shoulders compatible with S-Type masks

WingFrame L is a head and head-shoulder immobilisation system designed for use in radiotherapy treatments. It is compatible with S-type head and head-shoulders masks. The device can optionally include a shoulders retractor with 10 positions to adapt to the patient's needs.

Main features:

- Very lightweight and made of composite material for minimal attenuation.
- Compatible with standard TIMO (Foam) and SILVERMAN (PVC) headrests.
- Compatible with standard S-type head, head-neck, and head-shoulders masks.
- Optional shoulders retractor with 10 positions.
- Available in MR CONDITIONAL version up to 3 Tesla (WingFrame LfV-MR).

Product	Reference
Immobilisation System	
WingFrame LfC Carbon Fiber	AG1001076
WingFrame LfV-MR (MR compatible)	AG1001054
Accessories	
Shoulders retractor	AG1001073



*WingFrame LfC with shoulders retractor.